MOPEDIA

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I Class Medical Device according
to Directive 93/42/EEC and further modifications

1. Introduction

Thank you for purchasing a MOPEDIA swivel bathchair by Moretti S.p.A. This user's manual provides some suggestions as how to correctly use the product you have chosen and gives some valuable advice for your safety. Please read through the manual carefully before using the product. Should you have any queries, please contact your retailer for advice and assistance.

2. Intended use

Bath chair aids are intended to help people with walking difficulties to sit in the bath

3. EC Declaration of conformity

Moretti S.p.A. hereby declares, under its own exclusive responsibility, that the products manufactured and sold by Moretti S.p.A. and belonging to the cane aids category are in compliance with the applicable provisions outlined in the Medical Devices Directive 93/42/EEC and further modifications. To this end Moretti S.p.A. guarantees and declares as follows, under its own responsibility:

1. The devices in question satisfy the essential requirements of Annex I of directive 93/42/ EEC as prescribed by Annex VII of the above directive
2. The devices in question **ARE NOT MEASUREMENT DEVICES**
3. The devices in question **ARE NOT INTENDED FOR CLINICAL INVESTIGATION**
4. The devices in question are sold in **NON STERILE** packing
5. The devices in question should be considered as belonging to class I
6. **DO NOT** use/install the devices in question for any other purpose than that expressly indicated by Moretti S.p.A.
7. Moretti S.p.A. maintains and submits the documentation certifying conformity with Directive 93/42/EEC to the proper authorities and further modifications

4. General warnings

- Read carefully the entire manual before using the device
- Make sure that the product is delivered complete and keep it packed away from heat sources
- Service life- the device use limit is defined by the wear of the parts
- Do not allow children to play or operate the device

5. General description

5.1 Bathchair aids

RS938 Turning bathchair

The turning bathchair is manufactured from 1" steel which is powder coated to prevent rust. The seat is injection moulded plastic seat.

6. Before use

- Check if the frame is damaged in order to guarantee a safe use for people and things
- Check the correct assembly of the device

7. Assembly

7.1 Fitting the base frame over the bath

Rest the base frame securely on the bath. At least 30mm of each arm must rest on the sides of the bath. Turn the positioning screws on both sides until the black rubber ends are within 6mm

of the sides of the bath but not touching the sides. See Fig.1. this prevents the base frame from slipping.

7.2 Fitting the seat

Place the seat on the base frame. The socket under the seat locates over the swivel pin of the base frame. Ensure that the seat rotates freely. Lock the seat by pushing the black knob on the locking lever gently down until the lock locates in the slots in the base frame. (See Fig.2 and 3)

8. Warnings for a correct use

- The usage of the seat is recommended for patient weighting 130 kg
- DO NOT use with patient who exceed the weight indicated in the product label and in the technical features chapter of the present manual
- Avoid to leave the device in contact with water for long periods of time
- Periodically check the stability of the product, the screws are properly tightened and the rubber protection are in good conditions
- In case of doubts about use, adjustment or maintenance, contact your dealer
- Keep it away from children to avoid any possible risks
- The manufacturer and distributor decline all responsibility for any consequences resulting from an incorrect use of this product
- Check that the product has been delivered in perfect condition before you use it
- Worn, damaged or missing components must be replaced
- Help people with movement difficulties to seat on the bathchair
- Avoid any hard movement during patient positioning on the chair
- Pay attention on wet floors, can cause problems during patient positioning on the chair

9. How to use

- Turn the bathchair in front of the patient
- Help patient to take correct place on the seat
- Turn again the bathchair to put the patient in the bathtub

10. Maintenance

The MOPEDIA devices by Moretti are checked carefully and supplied with EC mark, once launched on the market. For patient and physician's safety we recommend to let check the device by the manufacturer or by laboratory authorized every 2 years. In case of repair, use only original spare parts and accessories.

11. Cleaning and disinfection

11.1 Cleaning

Use a damp cloth and a normal detergent. We recommend not to submerge the seat in the bath water, drain the seat tubes after use and dry the seat thoroughly. Don't use harsh abrasives or bleach based cleaners which may cause damage or discolouration.

11.2 Disinfection

If you need to disinfect the device use a common disinfectant cleaner

NOTE: Never use acids, alkalis or solvents such as acetone

12. Condition of disposal

12.1 General condition of disposal

In case of disposal of the device don't use the common container for municipal waste. We

recommend to dispose the device in the appropriate disposal areas for recycling.

13. Technical features

- Max. overall depth **50 cm**
- Max. overall width **72 cm**
- Max. overall height **54 cm**
- Seat width **42 cm**
- Seat depth **40 cm**
- Back support height **36 cm**
- Max weight **130 kg**

13. Warranty

Moretti products are guaranteed from material or manufacturing faults for a period of (two) 2 years from the purchasing date, subject to any limitations and exclusions specified below. The warranty does not apply in case of damages caused by improper use, abuse or modification of the product, and the warranty will not be valid if the instructions for use are not strictly followed. The correct intended use of the device is specified in this manual. Moretti is not responsible for damages resulting from personal injury or otherwise caused by or in relation to wrong installation or improper use. Moretti Moretti does not guarantee the products against damage or defects in the following conditions: natural disaster, not authorized maintenance or repair , damage due to power supply problems (where provided) use of spare parts not supplied by Moretti, improper use, not authorized alteration, shipping damage (other than original shipment by Moretti), or in case of insufficient maintenance as indicated in the manual. The warranty doesn't cover components subject to wear and tear during the correct use of the device.

13.1 Repairs

Warranty repairs In the event that a product Moretti to be defective in material or manufacturing faults during the warranty period, Moretti will confirm with customer if the fault can be covered from warranty. Moretti, in its sole discretion, can repair or replace the item under warranty at a Moretti dealer or at Moretti location . The labor cost for the repair may be charged to Moretti if the repair is covered by warranty. Repair or replacement does not renew or extend the warranty

Repair not covered by warranty A product not covered by warranty can be returned for repairing only after receiving prior authorization from Moretti customer service. The labor and shipping costs for a repair not covered by warranty will be fully borne by the customer or dealer. Repairs of products not covered by warranty are warranted for six (6) months from the date you receive the product repaired.

Non-defective device The customer will be informed if, after having examined and tested returned product, Moretti claims that the product is not defective. In this case the good will be returned to the customer, all shipping costs will be paid by customer.

13.2 Spare parts

The original Moretti spare parts are guaranteed for 6 (six) months from the receiving date.

13.3 Exempt clauses

Except as expressly specified in this warranty and within the limits of the law, Moretti does not offer any representation, warranty or condition, express or implied, including any

representations, warranties or conditions of merchantability, fitness for a particular purpose, non-infringement and interference. Moretti does not warrant that the use of the product Moretti will be uninterrupted or error free. The duration of any implied warranties that may be imposed by law are limited to the warranty period, within the limits of the law. Some states do not allow limitations on how long an implied warranty lasts or the exclusion or limitation of incidental or consequential damages in connection with products for consumers. In such states or countries, certain disclaimers and limitations of this warranty may not apply to you. This warranty is subject to change without notice.

WARRANTY CERTIFICATE

Product _____

Purchased on (date) _____

Retailer _____

Address _____ **Town/city** _____

Sold to _____

Address _____ **Town/city** _____

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11
www.morettispa.com email: info@morettispa.com
MADE IN TAIWAN

Note

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Via Bruxelles, 3 - Meleto
52022 Cavriglia (Arezzo)
Tel. +39 055 96 21 11

Fax. +39 055 96 21 200
www.morettispa.com
info@morettispa.com